



User's Manual

# Tissue Type Plasminogen Activator Human ELISA Kit



DEIA1625



96T



This product is for research use only and is not intended for diagnostic use.


For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For the quantitative measurement of Human Tissue type Plasminogen Activator concentrations in serum, plasma, tissue extracts, cell culture supernatants, milk, saliva and urine

### Storage

All kit components of this kit are stable at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

0.031 ng/mL-2 ng/mL

### Sensitivity

0.03 ng/mL