



**User's Manual**

# **ANA-8-Profile Antibody ELISA Kit**



**DEIA1678**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The ANA Combi assay is an indirect solid phase enzyme immunometric assay (ELISA) for the qualitative determination of ANA autoantibodies directed against RNP-70, RNP/Sm, Sm, SS-A, SS-B, Scl-70, Centromere B and Jo-1 in human serum or plasma.

### Storage

Store between 2-8°C. For more detailed information, please download the following document on our website.