



## User's Manual

# ssDNA Antibody ELISA Kit



DEIA1682



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

ssDNA Antibody ELISA kit is a solid phase enzyme immunoassay (ELISA) for the quantitative measurement of IgG class antibodies against single-stranded DNA in human serum or plasma.

### Storage

Store between 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

0-200 U/mL

