



User's Manual

Centromer B Antibody IgG ELISA Kit

REF

DEIA1684



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Anti-Centromer B is an indirect solid phase enzyme immunometric assay (ELISA). It is designed for measurement of IgG class autoantibodies directed against Centromer B.

Storage

Store between 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0-300 U/mL