



User's Manual

Estradiol ELISA Kit



DEIA1694



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Estradiol ELISA is an enzyme immunoassay for the quantitative in vitro diagnostic measurement of Estradiol in serum and plasma.

Storage

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. For more detailed information, please download the following document on our website.

Detection Range

9.7 pg/mL-2000 pg/mL.

Sensitivity

9.714 pg/mL

