



User's Manual

Human *Borrelia burgdorferi* IgG/IgM ELISA Kit



DEIA1717



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address:** 45-1 Ramsey Road, Shirley, NY 11967, USA

 **Tel:** 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  **Fax:** 1-631-938-8221

 **Email:** info@creative-diagnostics.com  **Web:** www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Borrelia burgdorferi IgG/IgM ELISA Test System is an enzyme-linked immunosorbent assay (ELISA) for the qualitative presumptive detection of total (IgG and IgM) antibodies to Borrelia burgdorferi in human serum.

Storage

Store components of the kit at 2-8°C or -20°C upon arrival up to the expiration date. For more detailed information, please download the following document on our website.