



User's Manual

Human Anti-Ovarian Antibody ELISA Kit



DEIA1725



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Anti-Ovarian Ab ELISA is a reliable test for measurement of antibodies directed against human ovarian tissues. This test is intended for the use with serum.

Storage

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. For more detailed information, please download the following document on our website.