



**User's Manual**

# Human Colony Stimulating Factor 2, CSF2 ELISA Kit



DEIA173



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This assay is a sandwich Enzyme Linked-Immuno-Sorbent Assay(ELISA). It is developed for quantitative measurement of Human CSF2 in serum, plasma and other biological fluids.

### Storage

Store at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

31.25-2,000 pg/mL

### Detection Limit

7 pg/mL