



## User's Manual

# Swine IFN- $\gamma$ ELISA Kit



DEIA1737



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Swine Interferon- $\gamma$  (Sw IFN- $\gamma$ ) ELISA is to be used for the in vitro quantitative determination of Sw IFN- $\gamma$  in swine serum, plasma, buffered solution, or cell culture medium. The assay will recognize both natural and recombinant Sw IFN- $\gamma$ .

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0-500 pg/mL

### Sensitivity

< 2.0 pg/mL