



User's Manual

Mouse Interleukin-1beta ELISA Kit



DEIA1752



192T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Mouse Interleukin-1 β (Ms IL-1 β) ELISA is to be used for the quantitative determination of Ms IL-1 β in Mouse serum, plasma (EDTA), buffered solution, or cell culture medium. The assay will recognize both natural and recombinant Ms IL-1 β .

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

0-1,000 pg/mL

Sensitivity

< 7 pg/mL