



User's Manual

Mouse Insulin ultrasensitive ELISA Kit



DEIA1762



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Ultrasensitive Mouse Insulin ELISA provides a method for determination of insulin in mouse serum or plasma.

Storage

Store kit at 2-8°C or -20°C upon arrival up to the expiration date. Avoid multiple freeze/thaw cycles. For more detailed information, please download the following document on our website.

Precision

Each sample was analyzed in four replicates on 16 different occasions.

Detection Range

The detection limit is 0.025 µg/L

