



User's Manual

Respiratory Syncytial Virus IgM ELISA Kit



DEIA1790



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Respiratory Syncytial Virus (RSV) IgM Enzyme Immunoassay Kit provides materials for the measurement of IgM-class antibodies to Respiratory Syncytial Virus in serum.

Storage

Reagents must be stored at 2-8°C. Once the foil bag has been opened, care should be taken to close it tightly again. For more detailed information, please download the following document on our website.