



User's Manual

Nucleosomes Antibody ELISA Kit



DEIA1831



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The kit is an indirect solid phase enzyme immunoassay (ELISA) for measurement of IgG class autoantibodies against Nucleosomes in human serum or plasma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0-200 U/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com