



**User's Manual**

# Ovarian Antibody Ig- Classifying ELISA Kit



DEIA1848



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Ovarian Antibody ELISA Kit is a reliable test for measurement of immunoglobulin class specific antibodies directed against human ovaries. This test is intended for the use with serum.

### Storage

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. For more detailed information, please download the following document on our website.



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)