



User's Manual

Parainfluenza virus Antibody IgA ELISA Kit



DEIA1857



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Parainfluenza virus Antibody IgA ELISA Kit is intended for measurement of IgA class antibodies against Parainfluenza Virus IgA in human serum or plasma (citrate).

Storage

The reagents are stable up to the expiry date stated on the label when stored at 2-8°C. For more detailed information, please download the following document on our website.