



User's Manual

Human Tumor Necrosis Factor (TNF superfamily, member 2), TNF ELISA Kit



DEIA186



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This assay is a sandwich Enzyme Linked-Immuno-Sorbent Assay(ELISA). It is developed for quantitative measurement of Human TNF in serum, plasma and other biological fluids.

Storage

Store at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

31.25-4,000 pg/mL

Detection Limit

4 pg/mL