



**User's Manual**

# Transglutaminase Antibody IgG ELISA Kit



DEIA1861



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Transglutaminase Antibody IgG ELISA Kit is an indirect solid phase enzyme immunoassay (ELISA) for measurement of IgG class autoantibodies against tissue Transglutaminase (tTG) in human serum or plasma.

### Storage

Store the reagents at 2-8°C (36°F-46°F). For more detailed information, please download the following document on our website.

### Detection Range

0-200 U/mL



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)