



User's Manual

ENA-4-Profile Antibody IgG ELISA Kit



DEIA1873



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

ENA-4-Profile Antibody IgG ELISA is an indirect solid phase enzyme immunometric assay (ELISA). It is designed for determination of IgG class autoantibodies directed against the extractable nuclear antigens SS-A (Ro), SS-B (La), Sm and RNP/Sm in human serum or plasma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

