



User's Manual

Alpha-Fodrin Antibody IgG/IgA ELISA Kit



DEIA1874



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

Alpha-Fodrin Antibody IgG/IgA ELISA Kit is an indirect solid phase enzyme immunoassay (ELISA) for measurement of IgG and IgA class autoantibodies against alpha-Fodrin in human serum or plasma.

Storage

Store the kit at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0-100 U/mL

