



User's Manual

Rat / Mouse Insulin ELISA Kit



DEIA1897



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Rat/Mouse Insulin ELISA kit is used for the non-radioactive determination of insulin in mouse and rat sera.

Storage

Prior to use, all components in the kit can be stored up to 2 weeks at 2-8°C.

For longer storage (>2 weeks), freeze diluted Wash Buffer, Insulin Standards, Quality Controls, and Matrix Solution at

Minimize repeated freeze and thaw of the Insulin Standards, Quality Controls and Matrix Solution.

Unused strips should be resealed in the foil pouch with the desiccant provided and stored at 2-8°C.

Refer to expiration dates on all reagents prior to use.

Do not mix reagents from different kits unless they have the same lot numbers.