



User's Manual

IGF-1 ELISA Kit



DEIA1901



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The IGF-1 ELISA is an enzyme immunoassay for the quantitative in vitro diagnostic measurement of Insulin-like Growth Factor 1 (IGF-1) in serum.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

1.29-600 ng/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: info@creative-diagnostics.com