



User's Manual

Human Total Proinsulin ELISA Kit



DEIA1906



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This Human Total Proinsulin ELISA kit is used for the non-radioactive measurement of human total proinsulin in serum and plasma.

Storage

Prior to use, all components in the kit can be stored up to 2 weeks at 2-8°C. For more detailed information, please download the following document on our website.