



## User's Manual

# Human C-Peptide ELISA Kit

REF

DEIA1928



96T

RUO



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This C-peptide ELISA provides a method for the quantitative determination of human C-peptide in serum, plasma or urine.

### Storage

For more detailed information, please download the following document on our website.

### Detection Range

The detection limit is 15 pmol/L (0.045 µg/L) calculated as two standard deviations above the Calibrator 0.