



## User's Manual

# Anti-LC-1 IgG ELISA Kit

REF

DEIA1933



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The Anti-LC-1 IgG ELISA Kit is used for determination of IgG antibodies to formiminotransferase-cyclodeaminase in human serum or plasma for the detection of autoimmune hepatitis (AIH).

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)