



User's Manual

Human Cystatin C ELISA Kit

REF

DEIA1955



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Cystatin C (human) ELISA is a sandwich enzyme immunoassay for easurement of human cystatin C.

Storage

Store the complete kit at 2-8°C. Under these conditions, the kit is stable until the expiration date (see label on the box). For more detailed information, please download the following document on our website.