



User's Manual

Anti-Tissue Transglutaminase ELISA Kit



DEIA1970



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

Anti-Tissue-Transglutaminase is an indirect solid phase enzyme immunoassay (ELISA) for the simultaneous measurement of IgG and IgA class autoantibodies to tissue transglutaminase (tTG) in human serum or plasma.

Storage

Store the kit reagents at 2-8°C. For more detailed information, please download the following document on our website.