



User's Manual

1,25-dihydroxy Vitamin D ELISA Kit

REF

DEIA1974



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The CD 1,25-Dihydroxy Vitamin D kit is a complete assay system intended for the purification of 1,25-dihydroxy Vitamin D (1,25D) in human serum or plasma by immunoextraction followed by quantitation by enzyme immunoassay.

Storage

Store the kit reagents at 2-8°C. For more detailed information, please download the following document on our website.

