



User's Manual

Human GIP ELISA Kit



DEIA1985



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Human GIP (Total) ELISA kit is used for the non-radioactive quantification of Human GIP in human serum, plasma, tissue extract and cell culture samples.

Storage

The kit should be stored at 2-8°C. For more detailed information, please download the following document on our website.