



## User's Manual

# Rat/Mouse GIP ELISA Kit



DEIA1986



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

This Rat/Mouse GIP (Total) ELISA kit is used for the non-radioactive quantification of Rat/Mouse GIP in rat or mouse serum, plasma, tissue extract and cell culture samples.

### Storage

Recommended storage for kit components is 2-8°C.

All components are shipped and stored at 2-8°C. Reconstituted standards and controls can be frozen for future use but repeated freeze thaws should be avoided.

Refer to expiration dates on all reagents prior to use.

Do not mix reagents from different kits unless they have the same lot numbers.

### Precision

#### Intra-assay Variation

Sample No.	Mean GIP Levels (pg/mL)	Intra-Assay %CV
1	92	1.0
2	98	2.4
3	173	5.9
4	150	1.1

#### Inter-assay Variation

Sample No.	Mean GIP Levels (pg/mL)	Inter-Assay %CV
1	80	5.9
2	113	1.1
3	178	4.1

### Sensitivity

8.2 pg/mL



Tel: 1-631-624-4882 (USA)

Tel: 44-161-818-6441 (Europe)



Fax: 1-631-938-8221



Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)