



User's Manual

Pancreatic Peptide ELISA Kit



DEIA2012



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This Enzyme Immunoassay kit is designed to measure a specific peptide and its related peptides based on the principle of "competitive" enzyme immunoassay.

Storage

Store the kit at 4°C upon receipt. For more detailed information, please download the following document on our website.

Detection Range

0.01-1000 ng/mL