



**User's Manual**

# Human HCG Beta free ELISA Kit



DEIA2027



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This product is an enzyme immunoassay for measurement of free  $\beta$ - subunit of human chorionic gonadotropin ( free  $\beta$ -hCG ) in serum.

### Storage

Store all contents at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

2.5-50 ng/mL

### Sensitivity

0.196 ng/mL