



User's Manual

Bovine Insulin ELISA Kit



DEIA2032



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This product provides a method for the quantitative determination of insulin in bovine serum and plasma.

Storage

The recommended storage temperature is 2-8°C.

Precision

Intra-assay Precision: 2.3-5.3%

Inter-assay Precision: 6.7-7.0%

Inter-Lot Precision: 6.8-7.3%

Detection Range

0.05-3 µg/L

Sensitivity

0.025 µg/L