



**User's Manual**

# **PTH-RP ( 1-34 ) ( human, rat, mouse ) ELISA Kit**



**DEIA2034**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

This Enzyme Immunoassay kit is designed to detect a specific peptide and its related peptides based on the principle of competitive enzyme immunoassay.

### Storage

Store the kit at 4°C upon receipt. Store 1X assay buffer at 4°C. Remove any unneeded strips from antibody-coated plate, reseal them in zip lock foil and keep at 4°C. Keep rehydrated solution of Standard, Biotinylated peptide, Antibody and HRP at 4°C.

### Detection Range

0-100 ng/mL