



User's Manual

Human sgp130 ELISA Kit



DEIA2035



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The Human sGP130 ELISA is to be used for measurement of gp130 protein (sGP130) in Human serum, plasma, buffered solutions or cell culture medium. The assay will recognize both natural and recombinant Human sGP130.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

56.25-1800 pg/mol