



**User's Manual**

# Human sRANKL (total) ELISA Kit



**DEIA2037**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Human sRANKL (total) ELISA is a sandwich enzyme immunoassay for measurement of total sRANKL (free and bound sRANKL) in serum and plasma samples.

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0.4-64 pmol/l

### Sensitivity

0.1 pmol/L