



User's Manual

Rat TNF-Alpha ELISA Kit



DEIA2038



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The rat TNF α ELISA is an enzyme-linked immunosorbent assay for the quantitative detection of rat TNF α .

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

39-2500 pg/mL

Sensitivity

11 pg/mL