



User's Manual

Human VEGF ELISA Kit



DEIA2045



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The VEGF ELISA is an enzyme-linked immunosorbent assay for measurement of Human VEGF levels (vascular endothelial growth factor) in cell culture supernatants, Human serum, plasma (EDTA, heparin and citrate) CSF and urine.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

9-700 pg/mL