



User's Manual

Mouse CRP ELISA Kit



DEIA2048



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The product is intended for the detection and quantification of mouse C-reactive protein (CRP) in mouse serum, plasma, and fluids.

Storage

The components should be stored at 2-8°C. Bring them to room temperature (20-25°C) before opening bottles and plate pouches.

Precision

Intra-assay Precision: 3-6%

Inter-assay Precision: 5-8%

Detection Range

1.56-25 ng/mL

Sensitivity

2.5 ng/mL

