



## User's Manual

# Pig *Trichinella spiralis* Ab ELISA Kit



DEIA2055



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

For measurement of antibodies in swine sera, plasma, whole blood or tissue fluid to *Trichinella spiralis* using the ELISA technique.

### Storage

Store between 2-8°C. Squeeze bottle containing diluted wash buffer may be stored at room temperature.