



User's Manual

Human sCD30 ELISA Kit

REF DEIA2057

Σ 96T



RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Human sCD30 ELISA Kit is an enzyme-linked immunosorbent assay for detection of Human sCD30.

Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

Detection Range

1.6-100 ng/mL

Sensitivity

0.33 ng/mL