



User's Manual

Human sCD40 ELISA Kit



DEIA2058



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

The human sCD40 ELISA Kit is an enzyme-linked immunosorbent assay for measurement of human sCD40.

Storage

Store kit reagents between 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

8-500 pg/mL

Sensitivity

7.92 pg/mL