



## User's Manual

# PAPP-A 60 ELISA Kit



DEIA2075



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

The PAPP-A 60 ELISA is an enzyme immunoassay for measurement of Pregnancy associated plasma protein A ( PAPP-A ) in serum and plasma.

### Storage

When stored at 2-8°C unopened reagents will retain reactivity until expiration date. Opened reagents must be stored at 2-8°C. Microtiter wells must be stored at 2-8°C. For more detailed information, please download the following document on our website.

### Detection Range

2.0-60.0 µg/mL

### Sensitivity

0.034 µg/mL

