



## User's Manual

# TriCat Research ELISA Kit



DEIA2077



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

### Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)  Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)

---

## PRODUCT INFORMATION

### Intended Use

Enzyme Immunoassay for measurement of Adrenaline ( Epinephrine ), Noradrenaline ( Norepinephrine ) and Dopamine. Flexible test system for various biological sample types and volumes.

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0.45-45 ng/mL

### Sensitivity

0.9-90 ng/mL