



User's Manual

Mouse Osteopontin ELISA Kit

REF

DEIA2089



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

This kit is to be used for the in-vitro quantitative determination of Mouse Osteopontin (Mouse OPN) in plasma (EDTA), urine, or cell culture media.

Storage

Store all contents at 2-8°C.

Precision

Intra-assay Precision: 4.0-4.4%

Inter-assay Precision: 4.0-4.4%

Detection Range

1-64 ng/mL

Sensitivity

0.15 ng/mL