



User's Manual

Human Acetylcholine Receptor Auto antibody ELISA Kit



DEIA2091



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Acetylcholine Receptor Autoantibody (AChRAb) ELISA kit is intended for use by professional persons only, for the determination of AChRAb in human serum.

Storage

Store all contents at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

0.5-20 nmol/l

Sensitivity

0.23 nmol/L