



User's Manual

Serotonin FAST ELISA Kit



DEIA2144



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This product is enzyme immunoassay for the quantitative determination of Serotonin in serum, urine and platelets.

Storage

Store all reagents at 2-8°C. Acylation Concentrate: room temperature. For more detailed information, please download the following document on our website.

Detection Range

15-2500 ng/mL

Sensitivity

5 ng/mL

