



User's Manual

Human Anti- Gliadin Deaminated Peptide IgA ELISA Kit

REF

DEIA2152



96T

RUO

This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The kit is a solid phase enzyme immunometric assay (ELISA) designed for the measurement of IgA class antibodies directed against deamidated Gliadin peptides (DGP) in human serum or plasma.

Storage

Store all contents at 2-8°C. For more detailed information, please download the following document on our website.

Sensitivity

90.9 %