



User's Manual

Carbonyl Protein ELISA Kit



DEIA2175



96T



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 Address: 45-1 Ramsey Road, Shirley, NY 11967, USA

 Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)  Fax: 1-631-938-8221

 Email: info@creative-diagnostics.com  Web: www.creative-diagnostics.com

PRODUCT INFORMATION

Intended Use

This ELISA Kit is intended for measurement of protein carbonyls in biological samples such as EDTA-plasma, bronchoalveolar lavage fluid and cerebrospinal fluid, cell extracts and other soluble protein samples.

Storage

Store all contents at 2-8°C. For more detailed information, please download the following document on our website.

Detection Range

156-2500 pmol/mg

Sensitivity

20 pmol/mg

