



**User's Manual**

# **Rat/Mouse Osteoprotegerin ELISA Kit**



**DEIA2177**



**96T**



This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

---

**Creative Diagnostics**

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)**  **Web: [www.creative-diagnostics.com](http://www.creative-diagnostics.com)**

---

## PRODUCT INFORMATION

### Intended Use

The Kit is a sandwich ELISA intended for the quantitative determination of Osteoprotege serum, plasma, urine and cell culture supernatants. It is for research use only.

### Storage

Store the reagents at 2-8°C until expiration date. For more detailed information, please download the following document on our website.

### Detection Range

0-4000 pg/mL



**Tel:** 1-631-624-4882 (USA)

**Tel:** 44-161-818-6441 (Europe)



**Fax:** 1-631-938-8221



**Email:** [info@creative-diagnostics.com](mailto:info@creative-diagnostics.com)