



User's Manual

Prostaglandin E2 HS ELISA Kit



DEIA2194



96T





This product is for research use only and is not intended for diagnostic use.

For illustrative purposes only. To perform the assay the instructions for use provided with the kit have to be used.

Creative Diagnostics

 **Address: 45-1 Ramsey Road, Shirley, NY 11967, USA**

 **Tel: 1-631-624-4882 (USA) 44-161-818-6441 (Europe)**  **Fax: 1-631-938-8221**

 **Email: info@creative-diagnostics.com**  **Web: www.creative-diagnostics.com**

PRODUCT INFORMATION

Intended Use

The Prostaglandin E2 HS ELISA Kit is designed to measure RBP.

Storage

Store contents at -20°C. For more detailed information, please download the following document on our website.

Detection Range

12.5-400 pg/mL

Sensitivity

10.9 pg/mL